ES

# EIA CanAg CA242

REF

101-10

IVD

CE

Instrucciones de uso. 2009-06

EN	EXPLANATION OF SYMBOLS
BG	ОБЯСНЕНИЕ НА СИМВОЛИТЕ
CS	VÝZNAM SYMBOLŮ
DA	SYMBOLFORKLARING
DE	ERKLÄRUNG DER SYMBOLE
EL	ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ
ES	SIGNIFICADO DE LOS SÍMBOLOS
ET	SÜMBOLITE SELGITUS
FR	EXPLICATION DES SYMBOLES
HR	OBJAŠNENJE SIMBOLA
HU	JELMAGYARÁZAT
IT	SPIEGAZIONE DEI SIMBOLI
LT	SIMBOLIŲ PAAIŠKINIMAI
LV	SIMBOLU SKAIDROJUMS
NL	VERKLARING DER SYMBOLEN
NO	SYMBOLFORKLARING
PL	OBJAŚNIENIE SYMBOLI
PT	EXPLICAÇÃO DOS SÍMBOLOS
RO	SEMNIȚAȚIA SIMBOLURILOR
RU	ОБОЗНАЧЕНИЯ
SE	SYMBOLFÖRKLARING
SK	VÝZNAM SYMBOLOV
SL	RAZLAGA SIMBOLOV
SR	OBJAŠNENJE SIMBOLA
TR	SEMBOLLERİN AÇIKLAMALARI

Use By/Годно до/Použitelné do/  
Holdbar til/Verwendbar bis/  
Ημερομηνία λήξης/Fecha  
de caducidad/Kölblük kuni/  
Utiliser jusque/Rok valjanosti/  
Felhasználható/Utilizzare entro/  
Sunautoti iki/Izletot līdz/Houdbaar  
tot/Brukes innen/Użyç przed/  
Prazo de validade/Expirã la/  
Использовать до/Använd före/  
Použite né do/ Uporabno do/  
Upotrebljivo do/Son Kullanna Tarihi

LOT

Batch code/Номер на партида/  
Číslo šarže/Lotnummer/  
Chargenbezeichnung/Αριθμός  
Παρτίδας/Código de lote/Partii  
kood/Code du lot/Kod serije/  
Sarzsám/Codice del lotto/  
Partijos kodus/Partijas kods/Lot  
nummer/Partikode/Kod partii/  
Código do lote/Număr de lot/  
Номер лота/Lotnummer/Číslo  
šarže/Številka serije/Kod partije/  
Parti Kodu

Date of manufacture/Дата на производство/Datum výroby/  
Produktionsdato/Herstellungsdatum/  
Ημερομηνία παραγωγής/Fecha de fabricación/Valmistamise kuupäev/  
Date de fabrication/Datum proizvodnje/  
Gyártási idő/Data di produzione/  
Pagaminimo data/Ražošanas datums/  
Productiedatum/Fremstillingsdato/  
Data produkcji/Data de fabrico/Data fabricației/Дата производства/  
Tillverkningsdatum/Dátum výroby/Datum izdelave/Datum proizvodnje/Üretim tarihi

Temperature limitation/  
Температурни граници/  
Терлотни омеzeи/  
Temperaturbegrænsning/  
Temperaturbegrenzung/  
Περιορισμοί θερμοκρασίας/  
Limites de temperatura/  
Temperatuuri piirang/  
Limite de température/  
Temperaturno ograničenje/  
Hőmérsékletre vonatkozó korlátozás/  
Limiti di temperatura/  
Temperatūriniai apribojimai/  
Temperatūras ierobežojums/  
Temperaturbepæring/  
Temperaturbegrensninger/  
Temperaturey granicne/  
Limite de temperatura/  
Limite de temperatură/  
Температурный режим/  
Temperaturbegrænsning/  
Teplotné obmedzenie  
Omejitve temperature/  
Temperaturno ograničenje/  
Sıcaklık sınırlaması/

## IVD

In Vitro Diagnostic Medical Device/  
Медицински уред за диагностика  
ин витро/Лéкаřský přístroj pro diagnostiku in vitro/Medicinsk udstyr til in vitro-diagnostik/In-vitro-Diagnostikum/  
Ιατροτεχνολογικό προϊόν για διάγνωση  
In Vitro/Dispositivo médico para diagnóstico in vitro/In vitro diagnostiline meditsiineseade/Dispositif médical de diagnostic in vitro/Diagnostički medicinski uređaj In Vitro/In vitro orvosdiagnostikai eszköz/Dispositivo medico per test diagnostici in vitro/In Vitro Diagnostinė Medicinos Priemonė/  
Medicinska ierīce in vitro diagnostikai/  
In vitro-diagnostisch medisch instrument/  
In vitro diagnostisk medisinsk utstyr/  
Wyrób medyczny do diagnostyki in vitro/  
Dispositivo Médico de Diagnóstico In Vitro/Dispozitiv medical pentru diagnostic in vitro/Только для диагностики In Vitro/Endast för in vitro-diagnostik/  
Zdravotnička pomôcka na diagnostiku in vitro/In vitro diagnostični pripomoček/  
Diagnostički medicinski uređaj In Vitro/<96> testleri için yeterlilik içerir

Contains sufficient for <96> tests/Съдържа достатъчно количество за тестове <96>/Lze použít pro <96> testů/Ineholder tilstrækkeligt/Inhalt ausreichend für <96> Prüfungen/Πεξεχόμενο επαρκές για «96» εξετάσεις/Contenido suficiente para <96> ensayos/Kogusest piisab <96> testi läbiviimiseks/Contenu suffisant pour «96» tests/Sadržaj dovoljno za <96> testova/A doboz tartalma <96> vizsgálat elvégzéséhez elegendő/Contenuto sufficiente per «96» saggi/Turiny's skirtas atlikti <96> tyrimus/ Saturs pietiekams <96> testiem/Inhoud voldoende voor «96» testen/til «96» test/ Tilstrækkelig innhold for <96> prøver/  
Wystarczy na wykonanie <96> testów/  
Conteúdo suficiente para «96» ensaios/  
Conținut suficient pentru 96 de teste/  
Содержит достаточные количества для «96» определений/Innehåller tillräckligt till «96» antal tester/Obsah postačuje na tento počet testov: <96>/Vsebinsa zadostuje za <96> testov/Sadržina dovoljna za <96> testova/<96> testleri için yeterlilik içerir

## REF

Catalogue number/Каталожен номер/  
Katalogové číslo/Katalognummer/  
Bestellnummer/Αριθμός καταλόγου/  
Número de catálogo/Katalogi number/  
Numéro de catalogue/Kataloški broj/  
Katalógusszám/Numero di catalogo/  
Katalogo numeris/Numurs katalogā/  
Catalogusnummer/Katalognummer/  
Numer katalogowy/Número do catálogo/  
Număr de catalog/Номер по каталогу/  
Produktnummer/Katalógové číslo/  
Kataloška številka/Kataloški broj/  
Katalog numarası

Consult Instructions for Use/  
Прочетете инструкцията за  
употреба/Konzultujte s návodem  
k použití/Se brugsanvisning/Siehe  
Gebrauchsanweisung/Συμβουλευτείτε  
της Οδηγίες σχετικά με τη χρήση/  
Consulte las instrucciones de uso/  
Vt kasutusjuhendit/Consulter le mode  
d'emploi/Pročítajte upute za uporabu/  
Olvassa el a használati utasítást/  
Consultare le istruzioni per l'uso/Dél  
naudojimo žiūrėkite instrukcijas/Izlasiet  
lietošanas instrukciju/Raadpleeg de  
instructies voor gebruik/Les instruksene  
for bruk/Sprawdzić w instrukcji użycia/  
Consulte as Instruções de Utilização/  
Consultați instrucțiunile de utilizare/  
Обратитесь к инструкции по  
применению/Se bruksanvisning/  
Prečítajte si návod na používanie/  
Pročitajte uputstvo za upotrebu/  
Kullanım Talimatlarını Bakınız

Contents of kit/Съдържание на набора/  
Obsah sady/Kittets inhold/Inhalt des  
Kits/Περιεχόμενα του κιτ/Contenido  
del kit/Komplekt sisaldab/Contenu du  
kit/Sadržaj opreme/A készlet tartalma/  
Contenuto del kit/Rinkinio turinys/  
Komplekta saturs/Inhoud van de set/  
Settets innhold/Zawartość zestawu/  
Conteúdo do kit/Conținutul setului/  
Компоненты набора/Kit innehåll/  
Obsah súpravy/Vsebina kompleta/Sadržaj  
opreme/Kitin içindekiler

Biological risks/Биологическа  
опасност/Biológická rizika/Biologisk  
fare/Biologische Gefahren/Βιολογικοί  
κίνδυνοι/Riesgos biológicos/  
Bioloogilised ohud/Risques biologiques/  
Biolóskli rizici/Biológiai kockázatok/Rischi  
biologici/Biologinis pavojus/Biológiskais  
risks/Biologische risico's/Biologische  
risikoer/Zagroženie biologiczne/Riscos  
biológicos/ Biologisk risk/Pericole  
biologice/Биологическая опасность/  
Biologicky rizikové/Biologické riziká/  
Biolóskli rizici/Biyolojik riskler

Human/C човешки производ/Ľidské/  
Human/Human/ἄνθρωπος αναφοράς/  
Humano/Inimāritolu/Humaine/Ljudskog  
porjekla/Humán/Origine Umana/  
Žmogaus kilmės/Cilvēku izcelsmes/  
Human/Menneske/Ludzka/Humano/  
Origine umână/Человеческого  
происхождения/Human/Ludské/  
Humanega izvora/Ljudskog porekla/İnsan

From mouse/C миши производ/Myši/  
Fra mus/der Maus/από ποντίκι/de ratón/  
Hiirtelt/De souris/Mišijeg porjekla/  
Egérböli/Murino/Pelės kilmės/No peles/  
Van muizen/Fra mus/Mysia/Do rato/De  
la șoareci/Мышиного происхождения/  
Från mus/Myšije/Mišijega izvora/Mišijeg  
porekla/Fareden

Bovine/C говежди производ/  
Hovēži/Bovin/Rind/από βοοειδή/  
Bovino/Veistelt/Bovine/Rogate stoke/  
Szarvasmarha/Bovino/Jaučio/No  
liellopa/Bovien/Bovini/Wolowy/Bovino/  
Origine bovină/крупного рогатого  
скота/Från ko/Hovädzie/Rogaveja  
izvora/Rogate krupne stoke/Bovin

Reconstitute with/Разтворение с/  
Rozfeďte pomocí/Rekonstitueres med/  
Rekonstituieren mit/Ανασύσταση με/  
Reconstituir con/Lahjendamine/  
Rekonstituer avec/Rekonstituiraite s/  
Feloldáshoz/Ricostituire con/Atkurti,  
ištirpdant su/Atšķaidīt ar/Rekonstituie  
met/Rekonstitueres med/Odtworzyć  
za pomocą/Reconstituir com/A  
se reconstitui cu/Разтворить в/  
Rekonstituera med/Rozriedte pomocou/  
Rekonstituiraite z/s/Ponovno formiranje  
sa/Yeniden oluşturulur

Manufacturer/Производитель/Výrobce/  
Producent/Hersteller/Κασκευαστής/  
Fabricante/Tootja/Fabricant/Proizvođač/  
Gyártó/Fabbricante/Gamintojas/  
Ražotājs/Fabrikant/Producent/  
Producent/Fabricante/Producător/  
Производитель/Тilverkare/ Výrobca/  
Izdelovalec/Proizvođač/Üretici

# EIA CanAg CA242

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1. Johansson, C., Nilsson, O., Bäckström, D., Jansson, E.-L. and Lindholm, L. (1991) Novel Epitopes on the CA50-Carrying Antigen: Chemical and Immunochemical Studies, *Tumor Biol.*, 12, 159-179.
2. Johansson, C., Nilsson, O. and Lindholm, L., (1991) Comparison of Serological Expression of Different Epitopes on the CA50-Carrying Antigen CanAg, *Int. J. Cancer*, 48, 757-763.
3. Nilsson, O., Johansson, C., Glimelius, B., Persson B., Nørgaard-Pedersen, B., Andrén-Sandberg, Å., and Lindholm L. (1992): Sensitivity and specificity of CA242 in gastro-intestinal cancer. A comparison with CEA, CA50 and CA19-9. *Br J Cancer* 65, 215-221.
4. Dahlén U., Karlsson B., Lindholm L., Nilsson O., (1993) Development of an enzyme immuno-assay for determination of the tumour associated antigen CA242, *J. Tumor Marker Oncology* 8, 3, p 111.
5. Kawa, S., Tokoo, M., Hasebe, O., Hayashi, K., Imai, H., Oguchi, H., Kiyosawa, K., Furuta, S., and Homma, T.,(1994) Comparative study of CA242 and CA19-9 for the diagnosis of pancreatic cancer, *Br. J. Cancer*, 70, 481-486.
6. Von Kleist, Hesse, Kananeh, (1996) Comparative Evaluation of Four Tumor Markers, CA 242, CA 19-9, TPA, and CEA in Carcinomas of the Colon, *Anticancer Research* 16: 2325-2332.
7. Spila A, Ferroni P, Cosimelli M, D'Alessandro R, Abbolito MR, Mariotti S, Aloe S, Carone MD, Graziano F, Tedesco M, Martini F, Mancini R, Stigliano V, Roselli M, Guadagni F. (2001) Comparative analysis of CA 242 and CA 19-9 serum tumor markers in colorectal cancer patients. A longitudinal evaluation. *Anticancer Res.* Mar-Apr;21(2B):1263-70.
8. Carpelan-Holmström, M., Haglund, K., Lundin, J., Järvinen, H and Roberts, P., (1996) Pre-operative serum levels of CA 242 and CEA predict outcome in colorectal cancer. *Eur. J. Cancer* 32(7), 1156-1161.
9. Engarås B. (2003 ) Individual cutoff levels of carcinoembryonic antigen and CA 242 indicate recurrence of colorectal cancer with high sensitivity. *Dis Colon Rectum.* Mar;46(3):313-21.
10. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
11. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986).

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**CanAg®** es una marca registrada de Fujirebio Diagnostics AB
**Fujirebio Diagnostics AB**