ES

# EIA CanAg PSA

REF

340-10

IVD

CE<sub>0197</sub>

Instrucciones de uso. 2009-11

EN	EXPLANATION OF SYMBOLS
BG	ОБЯСНЕНИЕ НА СИМВОЛИТЕ
CS	VÝZNAM SYMBOLŮ
DA	SYMBOLFORKLARING
DE	ERKLÄRUNG DER SYMBOLE
EL	ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ
ES	SIGNIFICADO DE LOS SÍMBOLOS
ET	SÜMBOLITE SELGITUS
FR	EXPLICATION DES SYMBOLES
HR	OBJAŠNJENJE SIMBOLA
HU	JELMAGYARÁZAT
IT	SPIEGAZIONE DEI SIMBOLI
LT	SIMBOLIŲ PAAIŠKINIMAI
LV	SIMBOLU SKAIDROJUMS
NL	VERKLARING DER SYMBOLEN
NO	SYMBOLFORKLARING
PL	OBJAŚNIENIE SYMBOLI
PT	EXPLICAÇÃO DOS SÍMBOLOS
RO	SEMNIȚAȚIA SIMBOLURILOR
RU	ОБОЗНАЧЕНИЯ
SE	SYMBOLFÖRKLARING
SK	VÝZNAM SYMBOLOV
SL	RAZLAGA SIMBOLOV
SR	OBJAŠNJENJE SIMBOLA
TR	SEMBOLLERİN AÇIKLAMALARI

Use By/Годно до/Použitelné do/  
Holdbar til/Verwendbar bis/  
Ημερομηνία λήξης/Fecha  
de caducidad/Kölblik kuni/  
Utiliser jusque/Rok valjanosti/  
Felhasználható/Utilizzare entro/  
Sunautoti iki/Izletot līdz/Houdbaar  
tot/Brukes innen/Użyç przed/  
Prazo de validade/Expirã la/  
Использовать до/Använd före/  
Použite né do/ Uporabno do/  
Upotrebljivo do/Son Kullanna Tarihi

LOT

Batch code/Номер на партида/  
Číslo šarže/Lotnummer/  
Chargenbezeichnung/Αριθμός  
Παρτίδας/Código de lote/Partii  
kood/Code du lot/Kod serije/  
Sarzsám/Codice del lotto/  
Partijos kodus/Partijas kods/Lot  
nummer/Partikode/Kod partii/  
Código do lote/Număr de lot/  
Номер лота/Lotnummer/Číslo  
šarže/Številka serije/Kod partije/  
Parti Kodu

Date of manufacture/Дата на производство/Datum výroby/  
Produktionsdato/Herstellungsdatum/  
Hμερομηνία παραγωγής/Fecha de fabricación/Valmistamise kuupäev/  
Date de fabrication/Datum proizvodnje/  
Gyártási idő/Data di produzione/  
Pagaminimo data/Ražošanas datums/  
Productiedatum/Fremstillingsdato/  
Data produkcji/Data de fabrico/Data fabricației/Дата производства/  
Tillverkningsdatum/Dátum výroby/Datum izdelave/Datum proizvodnje/Üretim tarihi

Temperature limitation/  
Температурни граници/  
Теплотни омеzeи/  
Temperaturbegrænsning/  
Temperaturbegrenzung/  
Περιορισμοί θερμοκρασίας/  
Limites de temperatura/  
Temperatuuri piirang/  
Limite de température/  
Temperaturno ograničenje/  
Hőmérsékletre vonatkozó korlátozás/  
Limiti di temperatura/  
Temperatūriniai apribojimai/  
Temperatūras ierobežojums/  
Temperaturbepæring/  
Temperaturbegrensninger/  
Temperaturey granične/  
Limite de temperatură/  
Температурный режим/  
Temperaturbegrænsning/  
Teplotné obmedzenie  
Omejitve temperature/  
Temperaturno ograničenje/  
Sıcaklık sınırlaması/

## IVD

In Vitro Diagnostic Medical Device/  
Медицински уред за диагностика  
ин витро/Лéкаřský přístroj pro  
diagnostiku in vitro/Medicinsk udstyr til  
in vitro-diagnostik/In-vitro-Diagnostikum/  
Ιατροτεχνολογικό προϊόν για διάγνωση  
In Vitro/Dispositivo médico para  
diagnóstico in vitro/In vitro diagnostiline  
meditsiiniiseade/Dispositif médical  
de diagnostic in vitro/Diagnostički  
medicinski uređaj In Vitro/In vitro  
orvosdiagnosztikai eszköz/Dispositivo  
medico per test diagnostici in vitro/In  
Vitro Diagnostinė Medicinos Priemonė/  
Medicinska ierice in vitro diagnostikai/  
In vitro-diagnostisch medisch instrument/  
In vitro diagnostisk medisinsk utstyr/  
Wyrób medyczny do diagnostyki in vitro/  
Dispositivo Médico de Diagnóstico In  
Vitro/Dispozitiv medical pentru diagnostic  
in vitro/Только для диагностики In  
Vitro/Endast för in vitro-diagnostik/  
Zdravotnička pomôcka na diagnostiku in  
vitro/In vitro diagnostični pripomoček/  
Diagnostički medicinski uređaj In  
Vitro/<96> testleri için yeterlilik içerir

Contains sufficient for <96> tests/Съдържа  
достатъчно количество за тестове  
<96>/Lze použít pro <96> testů/Ineholder  
tilstrækkeligt/Inhalt ausreichend für <96>  
Prüfungen/Πεξεχόμενο επαρκές για  
«96» εξετάσεις/Contenido suficiente para  
<96> ensayos/Kogusest piisab <96> testi  
läbiviimiseks/Contenu suffisant pour "96"  
tests/Sadržji dovoljno za <96> testova/A  
doboz tartalma <96> vizsgálat elvégzéséhez  
elegendő/Contenuto sufficiente per "96"  
saggi/Turiny's skirtas atlikti <96> tyrimus/  
Saturis pietiekams <96> testiem/Inhoud  
voldoende voor "96" testen/til "96" test/  
Tilstrækkelig innhold for <96> prøver/  
Wystarczy na wykonanie <96> testów/  
Conteúdo suficiente para "96" ensaios/  
Conținut suficient pentru 96 de teste/  
Содержит достаточные количества для  
«96» определений/Innehåller tillräckligt  
till "96" antal tester/Obsah postačuje na  
tento počet testov: <96>/Vsebinsa zadostuje  
za <96> testov/Sadržina dovoljna za <96>  
testova/<96> testleri için yeterlilik içerir

## REF

Catalogue number/Каталожен номер/  
Katalogové číslo/Katalognummer/  
Bestellnummer/Αριθμός καταλόγου/  
Número de catálogo/Katalogi number/  
Numéro de catalogue/Kataloški broj/  
Katalógusszám/Numero di catalogo/  
Katalogo numeris/Numurs katalogā/  
Catalogusnummer/Katalognummer/  
Numer katalogowy/Número do catálogo/  
Număr de catalog/Номер по каталогу/  
Produktnummer/Katalógové číslo/  
Kataloška številka/Kataloški broj/  
Katalog numarası

Consult Instructions for Use/  
Прочетете инструкцията за  
употреба/Konzultujte s návodem  
k použití/Se brugsanvisning/Siehe  
Gebrauchsanweisung/Συμβουλευτείτε  
της Οδηγίες σχετικά με τη χρήση/  
Consulte las instrucciones de uso/  
Vt kasutusjuhendit/Consulter le mode  
d'emploi/Pročítajte upute za uporabu/  
Olvassa el a használati utasítást/  
Consultare le istruzioni per l'uso/Dél  
naudojimo žiūrėkite instrukcijas/Izlasiet  
lietošanas instrukciju/Raadpleeg de  
instructies voor gebruik/Les instruksene  
for bruk/Sprawdzić w instrukcji użycia/  
Consulte as Instruções de Utilização/  
Consultați instrucțiunile de utilizare/  
Обратитесь к инструкции по  
применению/Se bruksanvisning/  
Prečítajte si návod na používanie/  
Pročitajte uputstvo za upotrebu/  
Kullanım Talimatlarını Bakınız

Contents of kit/Съдържание на набора/  
Obsah sady/Kittets indhold/Inhalt des  
Kits/Περιεχόμενα του κιτ/Contenido  
del kit/Komplekt sisaldab/Contenu du  
kit/Sadržaj opreme/A készlet tartalma/  
Contenuto del kit/Rinkinio turinys/  
Komplekta saturs/Inhoud van de set/  
Settets innhold/Zawartość zestawu/  
Conteúdo do kit/Conținutul setului/  
Компоненты набора/Kit innehåll/  
Obsah súpravy/Vsebina kompleta/Sadržaj  
opreme/Kitin içindekiler

Biological risks/Биологическа  
опасност/Biológická rizika/Biologisk  
fare/Biologische Gefahren/Biolóγικoi  
kίνδυνοι/Riesgos biológicos/  
Biolooigilised ohud/Risques biologiques/  
Biolóškli rizici/Biológiai kockázatok/Rischi  
biologici/Biologinis pavojus/Biológiskais  
risks/Biologische risico's/Biologiske  
risikoer/Zagroženie biologiczne/Riscos  
biológicos/ Biologisk risk/Pericole  
biologice/Биологическая опасность/  
Biologický rizikové/Biológické riziká/  
Biolóškli rizici/Biolojisk risiker

Human/C човешки производ/Lidské/  
Human/Human/δείγματα αναφοράς/  
Humano/Inimāritolu/Humaine/Ljudskog  
porjekla/Humán/Origine Umana/  
Žmogaus kilmės/Cilvēku izcelsmes/  
Human/Menneske/Ludzka/Humano/  
Origine umână/Человеческого  
происхождения/Human/Ludské/  
Humanega izvora/Ljudskog porekla/İnsan

From mouse/C миши производ/Myši/  
Fra mus/der Maus/από ποντίκι/de ratón/  
Hiirtelt/De souris/Mišijeg porjekla/  
Egérböli/Murino/Pelės kilmės/No peles/  
Van muizen/Fra mus/Mysia/Do rato/De  
la șoareci/Мышиного происхождения/  
Från mus/Myšije/Mišjega izvora/Mišijeg  
porekla/Fareden

Bovine/C говежди производ/  
Hovēži/Bovin/Rind/από βοοειδή/  
Bovino/Veistelt/Bovine/Rogate stoke/  
Szarvasmarha/Bovina/Jaučio/No  
liellopa/Bovien/Bovina/Wolowy/Bovino/  
Origine bovină/крупного рогатого  
скота/Från ko/Hovädzie/Rogaveja  
izvora/Rogate krupne stoke/Bovin

Reconstitute with/Разтваряне с/  
Rozfeđe pomoći/Rekonstitues med/  
Rekonstituieren mit/Ανασύσταση με/  
Reconstituir con/Lahjendamine/  
Reconstituer avec/Rekonstituiraite s/  
Feloldáshoz/Ricostituire con/LT/Atškaidīt  
ar/Reconstitutie met/Rekonstituerees  
med/Odtworzyć za pomocą/Reconstituir  
com/A se reconstitui cu/Растворить в/  
Rekonstituera med/Rozriedte pomocou/  
Rekonstituiraite z/s/Поновно formiranje  
sa/Yeniden oluşturalur

Manufacturer/Производитель/Výrobce/  
Producent/Hersteller/Κατασκευαστής/  
Fabricante/Tootja/Fabricant/Proizvođač/  
Gyártó/Fabbricante/Gamintojas/  
Ražotājs/Fabrikant/Produsent/  
Producent/Fabricante/Producător/  
Производитель/Тилverkare/ Výrobca/  
Izdevalavec/Proizvođač/Üretici

# EIA CanAg PSA

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1. Wang MC, Valenzuela LA, Murphy GP, Chu TM (1979). Purification of a human prostate specific antigen. *Invest Urol* 17: 159–163.
2. Lilja H. (1985). A kallikrein-like serine protease in prostatic fluid cleaves the predominant seminal vesicle protein. *J Clin Invest* 76: 1899–1903.
3. Oesterling JE (1991). Prostate specific antigen: A critical assessment of the most useful tumor marker for adenocarcinoma of the prostate. *J Urology* 145: 907–923.
4. Lilja H., Christensson A., Dahlén U., Matikainen M-T, Nilsson O., Pettersson K., Lövgren T. (1991). Prostate-specific antigen in serum occurs predominantly in complex with  $\alpha_1$ -antichymotrypsin. *Clin Chem* 37: 1618–1625.
5. Christensson A., Björk T., Nilsson O., Dahlén U., Matikainen M-T., Cockett ATK, Abrahamsson PA, Lilja H. (1993). Serum prostate specific antigen complexed to  $\alpha_1$ -antichymotrypsin as an indicator of prostate cancer. *J Urology* 150: 100–105.
6. Oesterling JE., Cooner WH., Jacobsen SJ., Guess HA., Lieber MM. (1993). The influence of patient age on the serum prostate specific antigen concentration: An important clinical observation. *Urol Clin North Am* 20: 671–80, 1993a.
7. Nilsson O., Peter A. Andersson I., Nilsson K., Grundström B., and Karlsson B. (1997) Antigenic determinants of prostatespecific antigen (PSA) and development of assays specific for different forms of PSA. *Br J Cancer* 75(6): 789–797.
8. P Price C., Allard J., Davies G., Dawney A., J Duffy M., France M., Mandarino G., Milford Ward A., Patel B., Sibley P. and Sturgeon C. (2001) Pre-and post-analytical factors that may influence use of serum prostate specific antigen and its isoforms in a screening programme for prostate cancer. *Ann Clin Biochem* 38: 188–216.
9. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
10. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986).

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**Fujirebio Diagnostics AB**