CanAg CA125 EIA

REF 400-10

IVD

CE

Kasutusjuhend. 2010-07

EN	EXPLANATION OF SYMBOLS
BG	ОБЯСНЕНИЕ НА СИМВОЛИТЕ
CS	VÝZNAM SYMBOLŮ
DA	SYMBOLFORKLARING
DE	ERKLÄRUNG DER SYMBOLE
EL	ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ
ES	SIGNIFICADO DE LOS SÍMBOLOS
ET	SÜMBOLITE SELGITUS
FR	EXPLICATION DES SYMBOLES
HR	OBJAŠNJENJE SIMBOLA
HU	JELMAGYARÁZAT
IT	SPIEGAZIONE DEI SIMBOLI
LT	SIMBOLIŲ PAAIŠKINIMAI
LV	SIMBOLU SKAIDROJUMS
NL	VERKLARING DER SYMBOLEN
NO	SYMBOLFORKLARING
PL	OBJAŚNIENIE SYMBOLI
PT	EXPLICAÇÃO DOS SÍMBOLOS
RO	SEMNIIFICAȚIA SIMBOLURILOR
RU	ОБОЗНАЧЕНИЯ
SE	SYMBOLFÖRKLARING
SK	VÝZNAM SYMBOLOV
SL	RAZLAGA SIMBOLOV
SR	OBJAŠNJENJE SIMBOLA
TR	SEMBOLLERİN AÇIKLAMALARI

Use By/Годно до/Použitelné do/
Holdbar til/Verwendbar bis/
Ημερομηνία λήξης/Fecha
de caducidad/Kölblik kuni/
Utiliser jusque/Rok valjanosti/
Felhasználható/Utilizzare entro/
Sunaudoti iki/Izlietot līdz/Houdbaar
tot/Brukes innen/Użyty przed/
Prazo de validade/Expíra la/
Использовать до/Använd före/
Použite né do/ Uporabno do/
Upotrebljivo do/Son Kullanna Tarihi

LOT

Batch code/Номер на партида/
Číslo šarže/Lotnummer/
Chargenbezeichnung/Αριθμός
Παρτίδας/Código de lote/Partii
kood/Code du lot/Kod serije/
Sarzszzám/Codice del lotto/
Partijos kodas/Partijas kods/Lot
nummer/Partikode/Kod partii/
Código do lote/Număr de lot/
Номер лота/Lotnummer/Číslo
šarže/Številka serije/Kod partije/
Parti Kodu

Date of manufacture/Дата на производство/Datum výroby/Produktionsdato/Herstellungsdatum/Hμερομηνία παραγωγής/Fecha de fabricación/Valmistamise kuupäev/Date de fabrication/Datum proizvodnje/Gyártási idő/Data di produzione/Pagaminimo data/Ražošanas datums/Productiedatum/Fremstillingsdato/Data produkcji/Data de fabrico/Data fabricației/Дата производства/Tillverkningsdatum/Dátum výroby/Datum izdelave/Datum proizvodnje/Üretim tarihi

In Vitro Diagnostic Medical Device/Медицински уред за диагностика ин витро/Lékařský přístroj pro diagnostiku in vitro/Medicinsk udstyr til in vitro-diagnostik/In-vitro-Diagnostikum/Ιατροτεχνολογικό προϊόν για διάγνωση In Vitro/Dispositivo médico para diagnóstico in vitro/In vitro diagnostiline meditsiiniseade/Dispositif médical de diagnostic in vitro/Diagnostički medicinski uređaj In Vitro/In vitro orvosdiagnostikai eszköz/Dispositivo medico per test diagnostici in vitro/In Vitro Diagnostinė Medicinos Priemonė/Mediciniska ierīce in vitro diagnostikai/In vitro-diagnostisch medisch instrument/In vitro diagnostisk medisinsk utstyr/Wyrób medyczny do diagnostyki in vitro/Dispositivo Médico de Diagnóstico In Vitro/Dispozitiv medical pentru diagnostic in vitro/Только для диагностики In Vitro/Endast för in vitro-diagnostik/Zdravotnícka pomôcka na diagnostiku in vitro/In vitro diagnostični pripomoček/Diagnostički medicinski uređaj In Vitro/<96> testleri için yeterlilik içerir

Catalogue number/Каталожен номер/Katalogové číslo/Katalognummer/Bestellnummer/Αριθμός καταλόγου/Número de catálogo/Katalogi number/Numéro de catalogue/Kataloški broj/Katalógusszám/Numero di catalogo/Katalogo numeris/Numurs katalogā/Catalogusnummer/Katalognummer/Numer katalogowy/Número do catálogo/Număr de catalog/Номер по каталогу/Produktnummer/Katalógové číslo/Kataloška številka/Kataloški broj/Katalog numarası

Temperature limitation/Температурни граници/Teplotní omezení/Temperaturbegrænsning/Temperaturbegrenzung/Περιορισμοί θερμοκρασίας/Límites de temperatura/Temperatuuri piirang/Limite de temperatură/Temperaturno ograničenje/Hőmérsékletre vonatkozó korlátozás/Limiti di temperatura/Temperatūriniai apribojimai/Temperatūras ierobežojums/Temperatuurbepierking/Temperaturbegrensninger/Temperaturey graniczne/Limite de temperatura/Limite de temperatură/Температурный режим/Temperaturbegrænsning/Teplotné obmedzenie/Omejitvev temperature/Temperaturno ograničenje/Sıcaklık sınırlaması

Contains sufficient for <96> tests/Съдържа достатъчно количество за тестове <96>/Lze použít pro <96> testů/Inneholder tilstrækkeligt/Inhalt ausreichend für <96> Prüfungen/Περιεχόμενο επαρκές για «96» εξετάσεις/Contenido suficiente para <96> ensayos/Kogusest piisab <96> testi läbiviimiseks/Contenu suffisant pour «96» tests/Sadrží dovoljno za <96> testova/A doboz tartalma <96> vizsgálat elvégzéséhez elegendő/Contenuto sufficiente per «96» saggi/Turiny's skirtas atlikti <96> tyrimus/Satur's pietiekams <96> testiem/Inhoud voldoende voor «96» testen/til "96" test/Tilstrækkelig innhold for <96> prøver/Wystarczy na wykonanie <96> testów/Conteúdo suficiente para "96" ensaios/Conținut suficient pentru 96 de teste/Содержит достаточные количества для «96» определений/Innehåller tillräckligt till "96" antal tester/Obsah postačuje na tento počet testov: <96>/Vsebinsa zadostuje za <96> testov/Sadržina dovoljna za <96> testova/<96> testleri için yeterlilik içerir


Consult Instructions for Use/
 Прочетете инструкцията за
 употреба/Konzultujte s návodem
 k použití/Se brugsanvisning/Siehe
 Gebrauchsanweisung/Συμβουλευτείτε
 τις Οδηγίες σχετικά με τη χρήση/
 Consulte las instrucciones de uso/
 Vt kasutusjuhendit/Consulter le mode
 d'emploi/Pročitajte upute za uporabu/
 Olvassa el a használati utasítást/
 Consultare le istruzioni per l'uso/Dël
 naudojimo žiūrėkite instrukcijas/Izlasiet
 lietošanas instrukciju/Raadpleeg de
 instructies voor gebruik/Les instruksene
 før bruk/Sprawdzic w instrukcji użycia/
 Consulte as Instruções de Utilização/
 Consultați instrucțiunile de utilizare/
 Обратитесь к инструкции по
 применению/Se bruksanvisning/
 Prečítajte si návod na používanie/
 Pročitajte uputstvo za upotrebu/
 Kullanım Talimatlarını Bakınız


Contents of kit/Съдържание на набора/
 Obsah sady/Kittets indhold/Inhalt des
 Kits/Περιεχόμενα του κιτ/Contenido
 del kit/Komplekt sisaldab/Contenu du
 kit/Sadržaj opreme/A készlet tartalma/
 Contenido del kit/Rinkinio turinys/
 Komplekta saturs/Inhoud van de set/
 Settets innhold/Zawartość zestawu/
 Conteúdo do kit/Conținutul setului/
 Компоненты набора/Kit innehåll/
 Obsah súpravy/Vsebina kompleta/Sadržaj
 opreme/Kitin içindekiler


Biological risks/Биологическа
 опасност/Biologická rizika/Biologisk
 fare/Biologische Gefahren/Βιολογικοί
 κίνδυνοι/Riesgos biológicos/
 Bioloogilised ohud/Risques biologiques/
 Biološkli rizici/Biológiai kockázatok/Rischi
 biologici/Biologinis pavojus/Bioloģiskais
 risks/Biologische risico's/Biologiske
 risikoer/Zagrozenie biologiczne/Riscos
 biológicos/ Biologisk risk/Pericole
 biologice/Биологическая опасность/
 Biologicky rizikové/Biologické riziká/
 Biološkli rizici/Biyolojik riskler

Human/С човешки произход/Lidské/
 Humant/Human/δείγματα αναφοράς/
 Humano/Inimpäritolu/Humaine/Ljudskog
 porjekla/Humán/Origine Umana/
 Žmogaus kilmės/Cilvēku izcelsmes/
 Human/Menneske/Ludzka/Humano/
 Origine umană/Человеческого
 происхождения/Human/Ludské/
 Humanega izvora/Ljudskog porekla/Ínsan

From mouse/С миши произход/Мыši/
 Fra mus/der Maus/από ποντίκι/de ratón/
 Hiirtelt/De souris/Mišijeg porjekla/
 Egérból/Murino/Pelés kilmés/No peles/
 Van muizen/Fra mus/Mysia/Do rato/De
 la șoareci/Мышиного происхождения/
 Från mus/Myšie/Mišjega izvora/Mišijeg
 porekla/Fareden

Bovine/С говежди произход/
 Hovězí/Bovin/Rind/από βοοειδή/
 Bovino/Veistelt/Bovine/Rogate stoke/
 Szarvasmarha/Bovino/Jaučio/No
 liellopa/Bovien/Bovin/Wołowy/Bovino/
 Origine bovină/крупного рогатого
 скота/Från ko/Hovädzie/Govejega
 izvora/Rogate krupne stoke/Bovin

Reconstitute with/Разтваряне с/
 Rozřed'te pomocí/Rekonstitueres med/
 Rekonstituieren mit/Ανασύσταση με/
 Reconstituir con/Lahjendamine/
 Reconstituer avec/Rekonstituiraite s/
 Feloldáshoz/Ricostituire con/Atkurti,
 ištirpdant su/Atšķaidīt ar/Reconstitue
 met/Rekonstitueres med/Odtworzyć
 za pomocą/Reconstituir com/A
 se reconstitui cu/Разтворит в/
 Rekonstituera med/Rozried'te pomocou/
 Rekonstituiraite z/s/Ponovno formiranje
 sa/Yeniden oluřturulur

Manufacturer/Производител/Výrobce/
 Producent/Hersteller/Κτασκευαστής/
 Fabricante/Tootja/Fabricant/Proizvođač/
 Gyártó/Fabbricante/Gamintojas/
 Ražotājs/Fabrikant/Produsent/
 Producent/Fabricante/Producător/
 Производител/Тиллverkare/ Výrobca/
 Izdelovalec/Proizvođač/Üretici

CanAg CA125 EIA

Kasutusjuhend. 2010-07

Immunoensüümtest

96 määramiseks

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ALLIKAD

1. Bast R., Feeney M., Lazarus H., Nadler L., Colvin R., Knapp R., (1981). Reactivity of a monoclonal antibody with human ovarian carcinoma. *J Clin Invest* 68: 1331 - 1337.
2. O'Brien T, Raymond L, Bannon G, Ford D, Hardardottir H, Miller F, Quirk J (1991). New monoclonal antibodies identify the glycoprotein carrying the CA125 epitope. *Am J Obstet Gynecol* 165: 1857 - 1864.
3. Nilsson O, Jansson E-L, Dahlen U, Nilsson K, Nustad K, Hogberg T, Lindholm L (1994). In: *Current Tumor Diagnosis: Applications, Clinical Relevance, Research, Trends*. Ed R. Klapdor. Zuckswardt Verlag Munchen, pp 401-405.
4. Nustad K. et al., (1996) Specificity and affinity of 26 monoclonal antibodies against the CA125 antigen: First report from the ISOBM TD-1 Workshop. *Tumor Biol* 17:196-219.
5. Bonfrer J., Duffy M., Radtke, M., Segurado O., Torre G., van Dalen A., Zwirner M. (1999) Tumor markers in gynaecological cancers-EGTM recommendations. *Anticancer Research* 19:2785-2820.
6. Tuxen M., Soletormos, G., Rustin, G., Nelstrop A., Dombernowsky, P. (2000) Biological variation and analytical imprecision of CA125 in patients with ovarian cancer. *Scand J Clin Invest* 60:713-722.
7. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
8. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13 , August (1986).

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CanAg® on firma Fujirebio Diagnostics AB registreeritud kaubamärk Fujirebio Diagnostics AB
Fujirebio Diagnostics AB