# CanAg S100 EIA

REF 708-10

IVD

CE

Instructions d'utilisation. 2010-05

EN	EXPLANATION OF SYMBOLS
BG	ОБЯСНЕНИЕ НА СИМВОЛИТЕ
CS	VÝZNAM SYMBOLŮ
DA	SYMBOLFORKLARING
DE	ERKLÄRUNG DER SYMBOLE
EL	ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ
ES	SIGNIFICADO DE LOS SÍMBOLOS
ET	SÜMBOLITE SELGITUS
FR	EXPLICATION DES SYMBOLES
HR	OBJAŠNJENJE SIMBOLA
HU	JELMAGYARÁZAT
IT	SPIEGAZIONE DEI SIMBOLI
LT	SIMBOLIŲ PAAIŠKINIMAI
LV	SIMBOLU SKAIDROJUMS
NL	VERKLARING DER SYMBOLEN
NO	SYMBOLFORKLARING
PL	OBJAŚNIENIE SYMBOLI
PT	EXPLICAÇÃO DOS SÍMBOLOS
RO	SEMNIIFICAȚIA SIMBOLURILOR
RU	ОБОЗНАЧЕНИЯ
SE	SYMBOLFÖRKLARING
SK	VÝZNAM SYMBOLOV
SL	RAZLAGA SIMBOLOV
SR	OBJAŠNJENJE SIMBOLA
TR	SEMBOLLERİN AÇIKLAMALARI

Use By/Годно до/Použitelné do/  
Holdbar til/Verwendbar bis/  
Ημερομηνία λήξης/Fecha  
de caducidad/Kölblik kuni/  
Utiliser jusque/Rok valjanosti/  
Felhasználható/Utilizzare entro/  
Sunaudoti iki/Izlietot līdz/Houdbaar  
tot/Brukes innen/Użyty przed/  
Prazo de validade/Expiră la/  
Использовать до/Använd före/  
Použite 'né do/Uporabno do/  
Upotrebljivo do/Son Kullanna Tarihi

LOT

Batch code/Номер на партида/  
Číslo šarže/Lotnummer/  
Chargenbezeichnung/Αριθμός  
Παρτίδας/Código de lote/Partii  
kood/Code du lot/Kod serije/  
Sarzszzám/Codice del lotto/  
Partijos kodas/Partijas kods/Lot  
nummer/Partikode/Kod partii/  
Código do lote/Număr de lot/  
Номер лота/Lotnummer/Číslo  
šarže/Številka serije/Kod partije/  
Parti Kodu

Date of manufacture/Дата на производство/Datum výroby/  
Produktionsdato/Herstellungsdatum/  
Ημερομηνία παραγωγής/Fecha de fabricación/Valmistamise kuupäev/  
Date de fabrication/Datum proizvodnje/  
Gyártási idő/Data di produzione/  
Pagaminimo data/Ražošanas datums/  
Productiedatum/Fremstillingsdato/  
Data produkcyj/Data de fabrico/Data fabricației/Дата производства/  
Tillverkningsdatum/Dátum výroby/Datum izdelave/Datum proizvodnje/Üretim tarihi

Temperature limitation/  
Температурни граници/  
Teplotní omezení/  
Temperaturbegrænsning/  
Temperaturbegrenzung/  
Περιορισμοί θερμοκρασίας/  
Limites de temperatura/  
Temperatuuri piirang/  
Limite de température/  
Temperaturno ograničenje/  
Hőmérsékletre vonatkozó korlátozás/  
Limiti di temperatura/  
Temperatūriniai apribojimai/  
Temperatūras ierobežojums/  
Temperatuurbepërking/  
Temperaturbegrensninger/  
Temperatury graniczne/  
Limite de temperatura/  
Limite de temperatură/  
Температурный режим/  
Temperaturbegrænsning/  
Teplotné obmedzenie  
Omejitve temperature/  
Temperaturno ograničenje/  
Sıcaklık sınırlaması/

In Vitro Diagnostic Medical Device/  
Медицински уред за диагностика  
ин витро/Lékařský přístroj pro diagnostiku in vitro/Medicinsk udstyr til in vitro-diagnostik/In-vitro-Diagnostikum/  
Ιατροτεχνολογικό προϊόν για διάγνωση  
In Vitro/Dispositivo médico para diagnóstico in vitro/In vitro diagnostiline meditsiiniseade/Dispositif médical de diagnostic in vitro/Diagnostički medicinski uređaj In Vitro/In vitro orvosdiagnostikai eszköz/Dispositivo medico per test diagnostici in vitro/In Vitro Diagnostiné Medicinos Priemonė/Mediciniska ierīce in vitro diagnostikai/  
In vitro-diagnostisch medisch instrument/  
In vitro diagnostisk medisinsk utstyr/  
Wyrób medyczny do diagnostyki in vitro/  
Dispositivo Médico de Diagnóstico In Vitro/Dispozitiv medical pentru diagnostic in vitro/Только для диагностики In Vitro/Endast för in vitro-diagnostik/  
Zdravotnícka pomôcka na diagnostiku in vitro/In vitro diagnostični pripomoček/  
Diagnostički medicinski uređaj In Vitro/<96> testleri için yeterlilik içerir

Contains sufficient for <96> tests/Σδържа достатъчно количество за тестове <96>/Lze použít pro <96> testů/Ineholder tilstrækkeligt/Inhalt ausreichend für <96> Prüfungen/Περιεχόμενο επαρκές για «96» εξετάσεις/Contenido suficiente para <96> ensayos/Kogusest piisab <96> testi läbiviimiseks/Contenu suffisant pour «96» tests/Sadrži dovoljno za <96> testova/A doboz tartalma <96> vizsgálat elvégzéséhez elegendő/Contenuto sufficiente per «96» saggi/Turiny's skirtas atlikti <96> tyrimus/Saturs pietiekams <96> testiem/Inhoud voldoende voor «96» testen/til «96» test/Tilstrekkelig innhold for <96> prøver/Wystarczy na wykonanie <96> testów/Conteúdo suficiente para «96» ensaios/Conținut suficient pentru 96 de teste/Содержит достаточные количества для «96» определений/Innehåller tillräckligt till «96» antal tester/Obsah postačuje na tento počet testov: <96>/Vsebinsa zadostuje za <96> testov/Sadržina dovoljna za <96> testova/<96> testleri için yeterlilik içerir

Catalogue number/Каталожен номер/  
Katalogové číslo/Katalognummer/  
Bestellnummer/Αριθμός καταλόγου/  
Número de catálogo/Kataloogi number/  
Numéro de catalogue/Kataloški broj/  
Katalógusszám/Numero di catalogo/  
Katalogo numeris/Numurs katalogā/  
Catalogusnummer/Katalognummer/  
Numer katalogowy/Número do catálogo/  
Număr de catalog/Номер по каталогу/  
Produktnummer/Katalógové číslo/  
Kataloška številka/Kataloški broj/  
Katalog numarası

Consult Instructions for Use/  
Прочетете инструкцията за  
употреба/Konzultujte s návodem  
k použití/Se brugsanvisning/Siehe  
Gebrauchsanweisung/Συμβουλευτείτε  
τις Οδηγίες σχετικά με τη χρήση/  
Consulte las instrucciones de uso/  
Vt kasutusjuhendit/Consulter le mode  
d'emploi/Pročitajte upute za uporabu/  
Olvassa el a használati utasítást/  
Consultare le istruzioni per l'uso/Dél  
naudojimo žiūrėkite instrukcijas/Izlasiet  
lietošanas instrukciju/Raadpleeg de  
instructies voor gebruik/Les instruksene  
før bruk/Sprawdzić w instrukcji użycia/  
Consulte as Instruções de Utilização/  
Consultați instrucțiunile de utilizare/  
Обратитесь к инструкции по  
применению/Se bruksanvisning/  
Prečítajte si návod na používanie/  
Pročitajte uputstvo za upotrebu/  
Kullanım Talimatlarını Bakınız

## CONT

Contents of kit/Съдържание на набора/  
Obsah sady/Kittets indhold/Inhalt des  
Kits/Περιεχόμενα του κιτ/Contenido  
del kit/Komplekt sisaldab/Contenu du  
kit/Sadržaj opreme/A készlet tartalma/  
Contenuto del kit/Rinkinio turinys/  
Komplekta saturs/Inhoud van de set/  
Settets innhold/Zawartość zestawu/  
Conteúdo do kit/Conținutul setului/  
Компоненты набора/Kit innehåll/  
Obsah súpravy/Vsebina kompleta/Sadržaj  
opreme/Kitin içindekiler

Biological risks/Биологическа  
опасност/Biologická rizika/Biologisk  
fare/Biologische Gefahren/Βιολογικοί  
κίνδυνοι/Riesgos biológicos/  
Bioloogilised ohud/Risques biologiques/  
Biološkli rizici/Biológiai kockázatok/Rischi  
biologici/Biologinis pavojus/Bioloģiskais  
risks/Biologische risico's/Biologiske  
risikoer/Zagroženie biologiczne/Riscos  
biológicos/ Biologisk risk/Pericole  
biologice/Биологическая опасность/  
Biologicky rizikové/Biologické riziká/  
Biološkli rizici/Bijolojik riskler

ORIG HUM

Human/С човешки произход/Lidské/  
Human/Human/δείγματα αναφοράς/  
Humano/Inimpăritolu/Humaine/Ljudskog  
porjekla/Humán/Origine Umana/  
Žmogaus kilmės/Cilvēku izcelsmes/  
Human/Menneske/Ludzka/Humano/  
Origine umană/Человеческого  
происхождения/Human/Ludské/  
Humanega izvora/Ljudskog porekla/İnsan

ORIG MOU

From mouse/С миши произход/Мыши/  
Fra mus/der Maus/από ποντίκι/de ratón/  
Hiirtelt/De souris/Mišijeg porjekla/  
Egérból/Murino/Pelés kilmés/No peles/  
Van muizen/Fra mus/Mysia/Do rato/De  
la șoareci/Мышиного происхождения/  
Frán mus/Мыšie/Mišjega izvora/Mišijeg  
porekla/Fareden

ORIG BOV

Bovine/С говежди произход/  
Hovězí/Bovin/Rind/από βοειδή/  
Bovino/Veistelt/Bovine/Rogate stoke/  
Szarvasmarha/Bovino/Jaučio/No  
liellopa/Bovien/Bovin/Wołowy/Bovino/  
Origine bovină/крупного рогатого  
скота/Frán ko/Hovädzie/Govejega  
izvora/Rogate krupne stoke/Bovin

Reconstitute with/Разтваряне с/  
Rozřed'te pomocí/Rekonstitueres med/  
Rekonstituieren mit/Ανασύσταση με/  
Reconstituir con/Lahjendamine/  
Reconstituer avec/Rekonstituiraĵte s/  
Feloldáshoz/Ricostituire con/Atkurti,  
ištirpdant su/Atšķaidīt ar/Reconstitutie  
met/Rekonstitueres med/Odtworzyć  
za pomocą/Reconstituir com/A  
se reconstitui cu/Растворить в/  
Rekonstituera med/Rozried'te pomocou/  
Rekonstituiraĵte z/s/Ponovno formiranje  
sa/Yeniden oluşturulur

Manufacturer/Производитель/Výrobce/  
Producent/Hersteller/Κτασκευαστής/  
Fabricante/Tootja/Fabricant/Proizvođač/  
Gyártó/Fabbricante/Gamintojas/  
Ražotājs/Fabrikant/Produsent/  
Producent/Fabricante/Producător/  
Производитель/Тиллverkare/ Výrobca/  
Izdelovalec/Proizvođač/Üretici

# CanAg S100 EIA

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Kit de Test Immunoenzymatique

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## UTILISATION

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## REFERENCES BIBLIOGRAPHIQUES

1. Moore BW (1965) A soluble protein characteristic of the nervous system. *Biochem Biophys Res Commun* 19:739-744.
2. Zimmer DB et al., (1995) The S100 protein family history, function and expression. *Brain Res Bull* 37:417-429.
3. Heizmann CW et al., (2002) S100 proteins:structure, functions and pathology. *Front Biosci* 7:1356-1368.
4. Schäfer BW et al. (1995) Isolation of a YAC clone covering a cluster of nine S100 genes on human chromosome 1q21: rationale for a new nomenclature of the S100 calcium-binding protein family. *Genomics* 25:638-643.
5. Takahashi K. et al., (1984) Immunohistochemical study on the distribution of  $\alpha$  and  $\beta$  subunits of S-100 protein in human neoplasm and normal tissues. *Virchows Arch* 45:385-396.
6. Banfalvi T. et al., (2003) Use of serum S-100B protein levels to monitor the clinical course of malignant melanoma. *Eur J Cancer* 39:164-169.
7. Djureen-Mårtensson E., et al., (2001) Serum S-100b protein as a prognostic marker in malignant cutaneous melanoma. *J Clin Oncol* 19:824-831.
8. Hauschild A. et al., (1999) S100B protein detection in serum is significant prognostic factor in metastatic melanoma. *Oncology* 56:338-344.
9. Wunderlich MT., et al., (1999) Early Neurobehavioral Outcome after stroke is related to release of neurobiochemical markers of brain damage. *Stroke* 30:1190-1195.
10. Martens P et al., (1998) Serum S100 and neuron specific enolase for prediction of regaining consciousness after global cerebral ischemia. *Stroke* 29:2363-2366.
11. Rosén H. et al., (1998) Increased serum levels of the S-100 protein are associated with hypoxic brain damage after cardiac arrest. *Stroke* 29: 473-477.
12. Ingebrigtsen T. et al.,(2000) The clinical value of serum S-100 protein measurements in minor head injury: a Scandinavian multicentre study. *Brain Inj* 14:1047-1055
13. Michetti F and Gazzolo D. (2002) S100B protein in biological fluids: A tool for perinatal medicine *Clin Chem* 48:2097-2104
14. Anderson R. et al., (2001) High serum S100B levels for trauma patients without head injuries *Neurosurgery* 48: 1255-1258
15. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
16. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986).

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CanAg® est une Marque Déposée de Fujirebio Diagnostics AB
Fujirebio Diagnostics AB