CanAg S100 EIA

REF 708-10

IVD

CE

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EN	EXPLANATION OF SYMBOLS
BG	ОБЯСНЕНИЕ НА СИМВОЛИТЕ
CS	VÝZNAM SYMBOLŮ
DA	SYMBOLFORKLARING
DE	ERKLÄRUNG DER SYMBOLE
EL	ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ
ES	SIGNIFICADO DE LOS SÍMBOLOS
ET	SÜMBOLITE SELGITUS
FR	EXPLICATION DES SYMBOLES
HR	OBJAŠNJENJE SIMBOLA
HU	JELMAGYARÁZAT
IT	SPIEGAZIONE DEI SIMBOLI
LT	SIMBOLIŲ PAAIŠKINIMAI
LV	SIMBOLU SKAIDROJUMS
NL	VERKLARING DER SYMBOLEN
NO	SYMBOLFORKLARING
PL	OBJAŚNIENIE SYMBOLI
PT	EXPLICAÇÃO DOS SÍMBOLOS
RO	SEMNIIFICAȚIA SIMBOLURILOR
RU	ОБОНАЧЕНИЯ
SE	SYMBOLFÖRKLARING
SK	VÝZNAM SYMBOLOV
SL	RAZLAGA SIMBOLOV
SR	OBJAŠNJENJE SIMBOLA
TR	SEMBOLLERİN AÇIKLAMALARI

Use By/Годно до/Použitelné do/
Holdbar til/Verwendbar bis/
Ημερομηνία λήξης/Fecha
de caducidad/Kölblik kuni/
Utiliser jusque/Rok valjanosti/
Felhasználható/Utilizzare entro/
Sunaudoti iki/Izlietot līdz/Houdbaar
tot/Brukes innen/Użyć przed/
Prazo de validade/Expiră la/
Использовать до/Använd före/
Použite né do/ Uporabno do/
Upotrebljivo do/Son Kullanma Tarihi

LOT

Batch code/Номер на партида/
Číslo šarže/Lotnummer/
Chargenbezeichnung/Αριθμός
Παρτίδας/Código de lote/Partii
kood/Code du lot/Kod serije/
Sarzs szám/Codice del lotto/
Partijos kodas/Partijas kods/Lot
nummer/Partikode/Kod partii/
Código do lote/Număr de lot/
Номер лота/Lotnummer/Číslo
šarže/Številka serije/Kod partije/
Parti Kodu

Date of manufacture/Дата на производство/Datum výroby/
Produktionsdato/Herstellungsdatum/
Ημερομηνία παραγωγής/Fecha de fabricación/
Valmistamise kuupäev/
Date de fabrication/Datum proizvodnje/
Gyártási idő/Data di produzione/
Pagaminimo data/Ražošanas datums/
Productiedatum/Fremstillingsdato/
Data produkcyj/Data de fabrico/Data fabricației/Дата производства/
Tillverkningsdatum/Dátum výroby/Datum izdelave/Datum proizvodnje/Üretim tarihi

Temperature limitation/
Температурни граници/
Teplotní omezení/
Temperaturbegrænsning/
Temperaturbegrenzung/
Περιορισμοί θερμοκρασίας/
Limites de temperatura/
Temperatuuri piirang/
Limite de température/
Temperaturno ograničenje/
Hőmérsékletre vonatkozó korlátozás/
Limiti di temperatura/
Temperatūriņiai apribojimai/
Temperatūras ierobežojums/
Temperatuurbepërking/
Temperaturbegrensninger/
Temperatury graniczne/
Limite de temperatura/
Limite de temperatură/
Температурный режим/
Temperaturbegrænsning/
Teplotné obmedzenie
Omejitve temperature/
Temperaturno ograničenje/
Sıcaklık sınırlaması/

In Vitro Diagnostic Medical Device/
Медицински уред за диагностика
ин витро/Lékařský přístroj pro diagnostiku in vitro/Medicinsk udstyr til in vitro-diagnostik/
In-vitro-Diagnostikum/
Ιατροτεχνολογικό προϊόν για διάγνωση
In Vitro/Dispositivo médico para diagnóstico in vitro/
In vitro diagnostiline meditsiiniseade/Dispositif médical de diagnostic in vitro/
Diagnostički medicinski uređaj In Vitro/
In vitro orvosdiagnostikai eszköz/Dispositivo medico per test diagnostici in vitro/
In Vitro Diagnostinè Medicinos Priemonė/
Medicīniska ierīce in vitro diagnostikai/
In vitro-diagnostisch medisch instrument/
In vitro diagnostisk medisinsk utstyr/
Wyrób medyczny do diagnostyki in vitro/
Dispositivo Médico de Diagnóstico In Vitro/
Dispozitiv medical pentru diagnostic in vitro/
Только для диагностики In Vitro/
Endast för in vitro-diagnostik/
Zdravotnícka pomôcka na diagnostiku in vitro/
In vitro diagnostični pripomoček/
Diagnostički medicinski uređaj In Vitro/
<96> testleri için yeterlilik içerir

Contains sufficient for <96> tests/Σύдържа достатъчно количество за тестове <96>/Lze použít pro <96> testů/
Ineholder tilstrækkeligt/
Inhalt ausreichend für <96> Prüfungen/
Περιεχόμενο επαρκές για «96» εξετάσεις/
Contenido suficiente para <96> ensayos/
Kogusest piisab <96> testi läbiviimiseks/
Contenu suffisant pour «96» tests/
Sadrži dovoljno za <96> testova/
A doboz tartalma <96> vizsgálat elvégzéséhez elegendő/
Contenuto sufficiente per «96» saggi/
Turinys skirtas atlikti <96> tyrimus/
Saturis pietiekams <96> testiem/
Inhoud voldoende voor «96» testen/
til «96» test/
Tilstrækkelig innhold for <96> prøver/
Wystarczy na wykonanie <96> testów/
Conteúdo suficiente para «96» ensaios/
Conținut suficient pentru 96 de teste/
Содержит достаточные количества для «96» определений/
Innehåller tillräckligt till «96» antal tester/
Obsah postačuje na tento počet testov:
<96>/Vsebina zadostuje za <96> testov/
Sadržina dovoljna za <96> testova/
<96> testleri için yeterlilik içerir

Catalogue number/Каталожен номер/
Katalogové číslo/Katalognummer/
Bestellnummer/Αριθμός καταλόγου/
Número de catálogo/Kataloogi number/
Numéro de catalogue/Kataloški broj/
Katalógusszám/Numero di catalogo/
Katalogo numeris/Numurs katalogā/
Catalogusnummer/Katalognummer/
Numer katalogowy/Número do catálogo/
Număr de catalog/Номер по каталогу/
Produktnummer/Katalógové číslo/
Kataloška številka/Kataloški broj/
Katalog numarası

Consult Instructions for Use/
 Прочетете инструкцията за
 употреба/Konzultujte s návodem
 k použití/Se brugsanvisning/Siehe
 Gebrauchsanweisung/Συμβουλευτείτε
 τις Οδηγίες σχετικά με τη χρήση/
 Consulte las instrucciones de uso/
 Vt kasutusjuhendit/Consulter le mode
 d'emploi/Pročitajte upute za uporabu/
 Olvassa el a használati utasítást/
 Consultare le istruzioni per l'uso/Dél
 naudojimo žiūrėkite instrukcijas/Izlasiet
 lietošanas instrukciju/Raadpleeg de
 instructies voor gebruik/Les instruksene
 før bruk/Sprawdzić w instrukcji użycia/
 Consulte as Instruções de Utilização/
 Consultați instrucțiunile de utilizare/
 Обратитесь к инструкции по
 применению/Se bruksanvisning/
 Prečítajte si návod na používanie/
 Pročitajte uputstvo za upotrebu/
 Kullanım Talimatlarını Bakınız

CONT

Contents of kit/Съдържание на набора/
 Obsah sady/Kittets indhold/Inhalt des
 Kits/Περιεχόμενα του κιτ/Contenido
 del kit/Komplekt sisaldab/Contenu du
 kit/Sadržaj opreme/A készlet tartalma/
 Contenuto del kit/Rinkinio turinys/
 Komplekta saturs/Inhoud van de set/
 Settets innhold/Zawartość zestawu/
 Conteúdo do kit/Conținutul setului/
 Компоненты набора/Kit innehåll/
 Obsah súpravy/Vsebina kompleta/Sadržaj
 opreme/Kitin içindekiler

Biological risks/Биологическа
 опасност/Biologická rizika/Biologisk
 fare/Biologische Gefahren/Βιολογικοί
 κίνδυνοι/Riesgos biológicos/
 Bioloogilised ohud/Risques biologiques/
 Biološkli rizici/Biológiai kockázatok/Rischi
 biologici/Biologinis pavojus/Bioloģiskais
 risks/Biologische risico's/Biologiske
 risikoer/Zagroženie biologiczne/Riscos
 biológicos/ Biologisk risk/Pericole
 biologice/Биологическая опасность/
 Biologicky rizikové/Biologické riziká/
 Biološkli rizici/Bijoloogik riskler

ORIG HUM

Human/С човешки произход/Lidské/
 Humant/Human/δείγματα αναφοράς/
 Humano/Inimpăritolu/Humaine/Ljudskog
 porjekla/Humán/Origine Umana/
 Žmogaus kilmės/Cilvēku izcelsmes/
 Human/Menneske/Ludzka/Humano/
 Origine umană/Человеческого
 происхождения/Human/Ludské/
 Humanega izvora/Ljudskog porekla/İnsan

ORIG MOU

From mouse/С миши произход/Мыši/
 Fra mus/der Maus/από ποντίκι/de ratón/
 Hiirtelt/De souris/Mišijeg porjekla/
 Egérból/Murino/Pelés kilmés/No peles/
 Van muizen/Fra mus/Mysia/Do rato/De
 la șoareci/Мышиного происхождения/
 Från mus/Myšie/Mišjega izvora/Mišijeg
 porekla/Fareden

ORIG BOV

Bovine/С говяжди произход/
 Hovězí/Bovin/Rind/από βοειδή/
 Bovino/Veistelt/Bovine/Rogate stoke/
 Szarvasmarha/Bovino/Jaučio/No
 liellopa/Bovien/Bovin/Wołowy/Bovino/
 Origine bovină/крупного рогатого
 скота/Från ko/Hovädzie/Govejega
 izvora/Rogate krupne stoke/Bovin

Reconstitute with/Разтваряне с/
 Rozředte pomocí/Rekonstituertes med/
 Rekonstituieren mit/Ανασύσταση με/
 Reconstituir con/Lahjendamine/
 Reconstituer avec/Rekonstituiraite s/
 Feloldáshoz/Ricostituire con/Atkurti,
 ištirpdant su/Atšķaidīt ar/Reconstitutie
 met/Rekonstituertes med/Odtworzyć
 za pomocą/Reconstituir com/A
 se reconstitui cu/Растворить в/
 Rekonstituera med/Rozriedte pomocou/
 Rekonstituiraite z/s/Ponovno formiranje
 sa/Yeniden oluşturulur

Manufacturer/Производитель/Výrobce/
 Producent/Hersteller/Κτασκευαστής/
 Fabricante/Tootja/Fabricant/Proizvođač/
 Gyártó/Fabbricante/Gamintojas/
 Ražotājs/Fabrikant/Produsent/
 Producent/Fabricante/Producător/
 Производител/Тиллverkare/ Výrobca/
 Izdelovalec/Proizvođač/Üretici

CanAg S100 EIA

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96 tests

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LITERATURE REFERENCES

1. Moore BW (1965) A soluble protein characteristic of the nervous system. *Biochem Biophys Res Commun* 19:739-744.
2. Zimmer DB et al., (1995) The S100 protein family history, function and expression. *Brain Res Bull* 37:417-429.
3. Heizmann CW et al., (2002) S100 proteins:structure, functions and pathology. *Front Biosci* 7:1356-1368.
4. Schäfer BW et al. (1995) Isolation of a YAC clone covering a cluster of nine S100 genes on human chromosome 1q21: rationale for a new nomenclature of the S100 calcium-binding protein family. *Genomics* 25:638-643.
5. Takahashi K. et al., (1984) Immunohistochemical study on the distribution of α and β subunits of S-100 protein in human neoplasm and normal tissues. *Virchows Arch* 45:385-396.
6. Banfalvi T. et al., (2003) Use of serum S-100B protein levels to monitor the clinical course of malignant melanoma. *Eur J Cancer* 39:164-169.
7. Djureen-Mårtensson E., et al., (2001) Serum S-100b protein as a prognostic marker in malignant cutaneous melanoma. *J Clin Oncol* 19:824-831.
8. Hauschild A. et al., (1999) S100B protein detection in serum is significant prognostic factor in metastatic melanoma. *Oncology* 56:338-344.
9. Wunderlich MT., et al., (1999) Early Neurobehavioral Outcome after stroke is related to release of neurobiochemical markers of brain damage. *Stroke* 30:1190-1195.
10. Martens P et al., (1998) Serum S100 and neuron specific enolase for prediction of regaining consciousness after global cerebral ischemia. *Stroke* 29:2363-2366.
11. Rosén H. et al., (1998) Increased serum levels of the S-100 protein are associated with hypoxic brain damage after cardiac arrest. *Stroke* 29: 473-477.
12. Ingebrigtsen T. et al., (2000) The clinical value of serum S-100 protein measurements in minor head injury: a Scandinavian multicentre study. *Brain Inj* 14:1047-1055
13. Michetti F and Gazzolo D. (2002) S100B protein in biological fluids: A tool for perinatal medicine *Clin Chem* 48:2097-2104
14. Anderson R. et al., (2001) High serum S100B levels for trauma patients without head injuries *Neurosurgery* 48: 1255-1258
15. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
16. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986).

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CanAg[®] è un marchio registrato di Fujirebio Diagnostics AB
Fujirebio Diagnostics AB